



• ScintLogic MaxMix

Liquid Scintillation Cocktail - High Capacity

- High capacity
- Good efficiency
- Non toxic solvent
- Odourless
- Biodegradable
- No loss through polyethylene vials

Good counting efficiency for a wide range of samples. Sample uptake up to 4ml of sample in 10ml of cocktail. Direct counting of biological samples. Odourless - 5 litre container - No diffusion through plastic vials.



• Counting Efficiency

ml of sample/10ml of ScintLogic MaxMix	% ³ H Efficiency
nil	46 %
1.0	40 %
2.0	36 %
3.0	35 %
5.0	30 %
10.0	21 %

Designed for low to medium ionic strength solutions. Many biological solutions can be counted without pre-treatment.

Non-hazardous - no transport documentation required.

Order No.	Description	Qty
LABSC15020	ScintLogic MaxMix	10 Litres (2x5)

