



## • ScintLogic HiCount Liquid Scintillation Cocktail - High Efficiency

- Odourless
- High efficiency
- Non-toxic solvent
- Good capacity
- Biodegradable
- No loss through polyethylene vials

High counting efficiency for low to medium ionic strength solutions. Sample uptake up to 3ml of sample in 10ml of cocktail. Direct counting of biological samples. Odourless - 5 litre container - No diffusion through plastic vials.




## • Counting Efficiency

ml of water/10ml of ScintLogic HiCount	% <sup>3</sup> H Efficiency		
	Water	0.15M NaCl	0.1M NaOH
nil	49	49	49
1.0	45	45	44
2.0	43	43	43
4.0	-	36	-

Designed for low to medium ionic strength solutions. Many biological solutions can be counted without pre-treatment.

Non-hazardous - no transport documentation required.

Order No.	Description	Qty
LABSC10010	ScintLogic HiCount	10 Litres (2x5)



Legend:  clear solution     two phase unsuitable for counting